

tag dental



INNOVATION.  
DESIGN.  
PRODUCTION.



# 05-40 implants

Geometry & Surface	08
axis	13
massif	25
crestone	33
Packaging	38



# 41-58 superstructures


Healing Caps	42
Pre Restoration	46
Final Restoration	47
Overdenture	52

# 59-80 instruments

Surgical Kits	60
Prosthetic Kit	67
Holding Kits	68
Drills	70
Surgical Tools	74
Packaging	78




# 85-87 DigiTag



SyBone 82

# tag dental

T.A.G. Medical Corporation LTD Products is an innovator in the field of medical and dental equipment and implants. Idyllically located in Kibbutz Gaaton, among the rolling hills and avocado groves of Israel's Western Galilee region, T.A.G. Medical Products was established as a subsidiary of the metal machine workshop that was founded on the kibbutz in 1949. By 1992, the company had broadened development and production capabilities, and began blazing the trail towards its current position as a leader in the manufacturing of surgical instruments.

T.A.G. Medical Products established two professional departments:

## tag arthro

The Orthopedic Division: Produced sterile and non-sterile medical and surgical devices and implants for Arthroscopy, Orthopedic and Endoscopy.

## tag dental

The Dental Division: Offers high quality implantology products, including: Dental implants, Superstructures, instruments, and biological materials.




T.A.G. is ISO 13485:2003 certified including CAN/CSA & CMDCAS by LNE.

All products are carry the CE Mark according to Annex II (Notified Body – LNE)


The scope of certifications covers the design, development, manufacturing and marketing of Orthopaedic and Dental Implants, Disposable and Reusable Surgical Instruments, including Force Gauges for rehabilitation.

The TAG Dental implant system products are cleared for marketing in the USA\* (K143326)

Quality is everyone's responsibility at our company, and we invest time and effort to ensure that we consistently deliver high quality products



Our expert development team with its many years of dental experience pleased to introduce a new generation of implants based on the concept of platform switching:



All implants are made of biocompatible Titanium alloy Grade 23 (Ti 6AL 4V ELI) .

The stability of the implant is a critical factor that will determine the final osseointegration results. All implants are Titanium alloy designed for all bone types.


# AdvanTAGes

CONNECTION

ONE  
IMPLANT  
SYSTEM


Internal Hex.  
Same for All Implants

Ø3.30-6.0 mm






PLATFORM  
SWITCHING


REVERSE CROWN  
CONCAVE PROFILE




MACRO POROUS




axis




crestone





massif



tag

Zn	Mg	Na	AL	K	P	Cl	V	S	Si	Ca	N	O	Ti	C	area	Implant
		0.34	0.84				0.29	0.08		0.24	0.21	52.58	15.81	29.57	#1	
			0.63				0.54			0.31	0.84	53.54	16.32	27.86	#2	



axis

# axis



### CONNECTION

Internal hexagon 2.44mm  
 Preventing implant/abutment rotation  
 Provides positioning indication  
 Provide accurate and firm insertion using dedicated tools.

### TOP

4 Micro threads at the top of the implant decrease the stress in the crestal zone and increase the contact area in the cortical bone.

### GEOMETRY

Wide pitch progressive thread enables bone compression during insertion, ultimately improving the bone's volume support, and positively affecting cancellous bone.

### SURFACE

The micro surface achieved by blasting technology followed by acid etching


### APICAL

Self-tapping allows an exceptional cutting capabilities, and is responsible for increased stability during insertion as the implant carves, fills, and compresses the bone.



## axis IMPLANTS RANGE

3.3



IAI-33010 L-10mm  
 IAI-33011 L-11.5mm  
 IAI-33013 L-13mm  
 IAI-33016 L-16mm

3.75



IAI-37508 L-8mm  
 IAI-37510 L-10mm  
 IAI-37511 L-11.5mm  
 IAI-37513 L-13mm  
 IAI-37516 L-16mm

4.2



IAI-42006 L-6mm  
 IAI-37510 L-8mm  
 IAI-42010 L-10mm  
 IAI-42011 L-11.5mm  
 IAI-42013 L-13mm  
 IAI-42016 L-16mm

5.0





IAI-50006 L-6mm  
 IAI-50008 L-8mm  
 IAI-50010 L-10mm  
 IAI-50011 L-11.5mm  
 IAI-50013 L-13mm

6.0



IAI-60006 L-6mm  
 IAI-60008 L-8mm  
 IAI-60010 L-10mm  
 IAI-60011 L-11.5mm  
 IAI-60013 L-13mm





# 2.44<sub>MM</sub> ONE CONNECTION FITS ALL

We are proud to present our new Axis dia. 3.30mm that completes the TAG Implants line and finalizes the concept of same connection for all diameters. The narrow is a 2.44mm internal hex. Connection which is the same as all our other implants diameter connection.

Due to the special design of the implant:

- Average static load test showed an average of 666N
- Fatigue limit of the implants tested obtained was 340N with 5 million cycles.

The static and fatigue test performed on our 3.3mm implants at the Technion (Haifa Israel) demonstrated exceptional resistance to load and cycle compared to other implants companies.

Implant 3.3








Fig. 4. Static Loading Curve






3.3




IAI-33010 L-10 mm  
IAI-33011 L-11.5mm  
IAI-33013 L-13mm  
IAI-33016 L-16mm

Ø 3.3mm Implant Procedure




	Soft Bone		Dense Bone	
Drill Speed (RPM)	1200-1500	900-1200	700-900	500-700
Diameter	Ø 1.90	Ø 2.0	Ø 2.5	Ø 2.8

This procedure does not replace the professional judgment of the surgeon




3.75




IAI-37508 L-8mm  
IAI-37510 L-10mm  
IAI-37511 L-11.5mm  
IAI-37513 L-13mm  
IAI-37516 L-16mm

Ø 3.75mm Implant Procedure




	Soft Bone			Dense Bone	
Drill Speed (RPM)	1200-1500	900-1200	500-700	400-700	200-400
Diameter	Ø 1.90	Ø 2.0	Ø 2.8	C 3.2	Ø 3.75/4.2

This procedure does not replace the professional judgment of the surgeon




4.2




IAI-42006 L-6mm IAI-42008 L-8mm IAI-42010 L-10mm IAI-42011 L-11.5mm IAI-42013 L-13mm IAI-42016 L-16mm

Ø 4.2mm Implant Procedure




	Soft Bone				Dense Bone	
Drill Speed (RPM)	1200-1500	900-1200	500-700	400-700	400-600	200-400
Diameter	Ø1.90	Ø2.0	Ø2.8	C3.2	C3.8	Ø3.5/4.2

This procedure does not replace the professional judgment of the surgeon




5.0




IAI-50006 L-6mm IAI-50008 L-8mm IAI-50010 L-10mm IAI-50011 L-11.5mm IAI-50013 L-13mm

Ø 5.0mm Implant Procedure




	Soft Bone				Dense Bone	
Drill Speed (RPM)	1200-1500	900-1200	500-700	400-700	400-600	300-500
Diameter	Ø1.90	Ø2.0	Ø2.8	C3.2	C3.8	C4.5
						Ø5/6

This procedure does not replace the professional judgment of the surgeon



6.0



IAI-60006  
L-6mm


IAI-60008  
L-8mm

IAI-60010  
L-10mm

IAI-60011  
L-11.5mm

IAI-60013  
L-13mm

Ø 6.0mm Implant Procedure





Drill Speed (RPM)

Diameter

	Soft Bone			Dense Bone				
Drill Speed (RPM)	1200-1500	900-1200	500-700	400-700	400-600	300-500	450-250	200-400
Diameter	Ø1.90	Ø2.0	Ø2.8	C3.2	C3.8	C4.5	C5.5	Ø5/6

This procedure does not replace the professional judgment of the surgeon



MASSIF

# massif

One of our most advanced implants - a cylindrical screw-type implant



### CONNECTION

Internal hexagon 2.44mm  
Preventing implant/abutment rotation  
Provides positioning indication  
Provide accurate and firm insertion using dedicated tools.

### TOP

Micro rings - 4 micro rings and connecting micro thread at the top of the implant increase the envelop surface as well as bone to implant contact, thus improving the implant's primary stability and interfacial shear strength, especially in the crestal zone.

### GEOMETRY

A progressive thread allows bone compression and initial stability in the bone, and enables load transfer from the crestal to the cancellous bone.  
The double macro thread offers quick implantation, while the deep profile reverse buttress provides high bone grip and stability.

### SURFACE

The entire implant surface treatment includes Blasting technology and Acid Etching.

### APICAL

Self-tapping allows an exceptional cutting capabilities, and is responsible for increased stability during insertion as the implant carves, fills, and compresses the bone.




## massif IMPLANTS RANGE

3.75



4.2




5.0




6.0






3.75




IMI-37508 L-8mm  
 IMI-37510 L-10mm  
 IMI-37511 L-11.5mm  
 IMI-37513 L-13mm  
 IMI-37516 L-16mm

Ø 3.75mm Implant Procedure




	Soft Bone			Dense Bone	
Drill Speed (RPM)	1200-1500	900-1200	500-700	400-700	200-400
Diameter	Ø1.90	Ø2.0	Ø2.8	Ø3.2	Ø3.75/4.2

This procedure does not replace the professional judgment of the surgeon




4.2




IMI-42006 L-6mm  
 IMI-42008 L-8mm  
 IMI-42010 L-10mm  
 IMI-42011 L-11.5mm  
 IMI-42013 L-13mm  
 IMI-42016 L-16mm

Ø 4.2mm Implant Procedure




	Soft Bone				Dense Bone	
Drill Speed (RPM)	1200-1500	900-1200	500-700	400-700	400-600	200-400
Diameter	Ø1.90	Ø2.0	Ø2.8	Ø3.2	Ø3.8	Ø3.75/4.2


This procedure does not replace the professional judgment of the surgeon



5.0




Ø 5.0mm Implant Procedure




	Soft Bone					Dense Bone	
Drill Speed (RPM)	1200-1500	900-1200	500-700	400-700	400-600	300-500	200-400
Diameter	Ø1.90	Ø2.0	Ø2.8	Ø3.2	Ø3.8	Ø4.5	Ø5/6


This procedure does not replace the professional judgment of the surgeon



6.0





Ø 6.0mm Implant Procedure



	Soft Bone						Dense Bone		
Drill Speed (RPM)	1200-1500	900-1200	500-700	400-700	400-600	300-500	200-400	200-450	200-400
Diameter	Ø1.90	Ø2.0	Ø2.8	Ø3.2	Ø3.8	Ø4.5	Ø5	Ø5.5	Ø5/6

This procedure does not replace the professional judgment of the surgeon





crestone

# crestone

A one piece implant, designed with an integrated abutment for a single stage surgical procedure



### CONNECTION

External square 2.10mm  
Preventing implant/abutment rotation  
Provides positioning indication  
Provide accurate and firm insertion using dedicated tools.

### TOP

Crestone implant design features an integrated abutment for a single stage surgical procedure, and a force transferring geometry (2.1 square) for firm insertion.

### GEOMETRY

Innovative design with dual thread for easy insertion ensures maximum strength and stability.


### SURFACE

The micro surface achieved by blasting technology followed by acid etching.

### APICAL

Self-tapping allows for exceptional cutting capabilities, and is responsible for increased stability during insertion as the implant carves, fills, and compresses the bone.

Made of biocompatible Titanium Alloy Grade 23 (TI 6AL 4V ELI).  
It is recommended to use the Crestone in narrow ridges and tight spaces.  
Available in several Lengths and Diameters.




## crestone IMPLANTS RANGE

3.0



3.5




PP1-0002

One Piece Analog  
19mm / Stainless Steel

TK1-0011

Square Driver One Piece  
28mm / Stainless Steel




3.0




ICI-30010 L-10mm    ICI-30011 L-11.5mm    ICI-30013 L-13mm    ICI-30016 L-16mm

Ø 3.0mm Implant Procedure



	Soft Bone		Dense Bone	
Drill Speed (RPM)	1200-1500	900-1200	700-900	500-700
Diameter	Ø1.90	Ø2.0	Ø2.5	Ø2.8

This procedure does not replace the professional judgment of the surgeon




3.5



ICI-35010 L-10mm    ICI-35011 L-11.5mm    ICI-35013 L-13mm    ICI-35016 L-16mm

Ø 3.5mm Implant Procedure



	Soft Bone		Dense Bone	
Drill Speed (RPM)	1200-1500	900-1200	700-900	500-700
Diameter	Ø1.90	Ø2.0	Ø2.5	Ø2.8


This procedure does not replace the professional judgment of the surgeon

Drill Speed (RPM)  
Diameter



- Visual identification of implant diameter by cup color
- Implant length and diameter are labeled upon the cap
- Each implant is packed in a sealed sterile tube and closed within a shrink sleeve for absolute impermeability





# HEALING CAPS

## NARROW

3.75



PHI-0016  
Healing Cap  
2mm



PHI-0002  
Healing Cap  
3mm



PHI-0003  
Healing Cap  
4mm



PHI-0004  
Healing Cap  
5mm



PHI-0005  
Healing Cap  
6mm

## STANDARD

4.7



PHI-0017  
Anatomic  
Healing Cap  
2mm



PHI-0006  
Anatomic  
Healing Cap  
3mm



PHI-0007  
Anatomic  
Healing Cap  
4mm



PHI-0008  
Anatomic  
Healing Cap  
5mm



PHI-0009  
Anatomic  
Healing Cap  
6mm

## WIDE

5.5



PHI-0018  
Wide Anatomic  
Healing Cap 2mm



PHI-0013  
Wide Anatomic  
Healing Cap 3mm



PHI-0014  
Wide Anatomic  
Healing Cap 4mm



PHI-0015  
Wide Anatomic  
Healing Cap 5mm




PHI-0019  
Wide Anatomic  
Healing Cap 6mm

L6mm




Implant Abutment Analog  
PPI- 0001  
14mm / Stainless Steel

Abutment Analog


Implant Analog




# TRANSFERS



## Peek Esthetic Abutments



## Titanium Esthetic Abutments



0°



15°



25°



0°



15°



25°




## UNIVERSAL



All our abutments are supplied together with the relevant screw. The titanium Abutment should be tightened to 30-35 Ncm.



# Burn-Out Esthetic Abutments



0°

Direct Plastic

							
Direct Plastic Cylinder With Hexagon	Direct Plastic Cylinder Round	Plastic Cylinder with Hex. Ti	Plastic Cylinder Round Ti	Direct gold Plastic abutment Cylinder with hex.	Direct gold Plastic abutment Cylinder Round	Direct CoCr Plastic abutment Cylinder hex.	Direct CoCr Plastic Abutment Cylinder Round
L 11mm Acetal	L 11mm Acetal	L 15mm Acetal, Titanium	L 15mm Acetal, Titanium	L 15mm Acetal, Gold	L 15mm Acetal, Gold	L 15mm Acetal, CoCr	L 15mm Acetal, CoCr

0°

		
Esthetic Abutment H1mm	Esthetic Abutment H2mm	Esthetic Abutment H3mm
L 9mm / Acetal	L 9mm / Acetal	L 9mm / Acetal

15°

		
15° Esthetic Abutment H1mm	15° Esthetic Abutment H2mm	15° Esthetic Abutment H3mm
L 9mm / Acetal	L 9mm / Acetal	L 9mm / Acetal

25°


		
25° Esthetic Abutment H1mm	25° Esthetic Abutment H2mm	25° Esthetic Abutment H3mm
L 8mm / Acetal	L 8mm / Acetal	L 8mm / Acetal

# ANATOMIC RESTORATION KIT



4mm	PF1-ARK14 Anatomic Restoration Kit H1mm L4mm	PF1-ARK24 Anatomic Restoration Kit H2mm L4mm	PF1-ARK34 Anatomic Restoration Kit H3mm L4mm	PF1-ARK44 Anatomic Restoration Kit H4mm L4mm	PF1-ARK54 Anatomic Restoration Kit H5mm L4mm
6mm	PF1-ARK16 Anatomic Restoration Kit H1mm L6mm	PF1-ARK26 Anatomic Restoration Kit H2mm L6mm	PF1-ARK36 Anatomic Restoration Kit H3mm L6mm	PF1-ARK46 Anatomic Restoration Kit H4mm L6mm	PF1-ARK56 Anatomic Restoration Kit H5mm L6mm
8mm	PF1-ARK18 Anatomic Restoration Kit H1mm L8mm	PF1-ARK28 Anatomic Restoration Kit H2mm L8mm	PF1-ARK38 Anatomic Restoration Kit H3mm L8mm	PF1-ARK48 Anatomic Restoration Kit H4mm L8mm	PF1-ARK58 Anatomic Restoration Kit H5mm L8mm

# Anatomic Abutments



## 4mm

PF1-0005	PF1-0008	PF1-0011	PF1-0014	PF1-0087
Anatomic Abutment H1mm L 4mm Titanium	Anatomic Abutment H2mm L 4mm Titanium	Anatomic Abutment H3mm L 4mm Titanium	Anatomic Abutment H4mm L 4mm Titanium	Anatomic Abutment H5mm L4mm Titanium

## 6mm

PF1-0006	PF1-0009	PF1-0012	PF1-0015	PF1-0088
Anatomic Abutment H1mm L 6mm Titanium	Anatomic Abutment H2mm L 6mm Titanium	Anatomic Abutment H3mm L 6mm Titanium	Anatomic Abutment H4mm L 6mm Titanium	Anatomic Abutment H5mm L6mm Titanium

## 8mm

PF1-0007	PF1-0010	PF1-0013	PF1-0016	PF1-0089
Anatomic Abutment H1mm L 8mm Titanium	Anatomic Abutment H2mm L8mm Titanium	Anatomic Abutment H3mm L 8mm Titanium	Anatomic Abutment H4mm L 8mm Titanium	Anatomic Abutment H5mm L8mm Titanium

## 6mm Wide Platform


PF1-0065	PF1-0064	PF1-0063	PF1-0090
Anatomic Abutment H2mm L 6mm Wide Titanium	Anatomic Abutment H3mm L 6mm Wide Titanium	Anatomic Abutment H4mm L 6mm Wide Titanium	Anatomic Abutment H5mm L6mm Wide Titanium




All our abutments are supplied together with the relevant screw. The titanium Abutment should be tightened to 30-35 Ncm.

# OVERDENTURE

 PATENT PENDING



**Ball Attachment Angulated ( BAA) - Patent Pending**  
The BAA system provides the tools to restore even compromised edentulous cases.  
Made in one piece the BAA is an esthetic design.  
The Reduced size will improve the esthetic restoration.



- 1 One-piece Design
- 2 Reduced Size Design
- 3 Intelligent Design which Allows for Restorative Flexibility



**Multi Unit Angulated (MUA) - Patent Pending**  
Made in one piece the MUA have is an esthetic and compact design.  
Reduced size improves the esthetic restoration.

MUA and BAA are intended for use with implants of diameter 3.75mm -6.0mm

# MULTI-UNIT ABUTMENTS



POI-0016	POI-0017	POI-0018	POI-0019	POI-0020		
Multi-Unit H 1mm	Multi-Unit H 2mm	Multi-Unit H 3mm	Multi-Unit H 4mm	Multi-Unit H 5mm		
Titanium	Titanium	Titanium	Titanium	Titanium		
POI-0038	POI-0039	POI-0047	POI-0040	POI-0041	POI-0049	TKI-0033
Multi-Unit 17° / H 2mm	Multi-Unit 17° / H 3mm	Multi-Unit 17° / H 4mm	Multi-Unit 30° / H 2mm	Multi-Unit 30° / H 3mm	Multi-Unit 30° / H 4mm	Key for M.U&balls
Titanium	Titanium	Titanium	Titanium	Titanium	Titanium	Stainless Steel

PATENT PENDING


# MULTI - UNIT ACCESSORIES

POI-0027	POI-0028				
Multi-Unit Plastic Cylinder	Multi-Unit Temporary Ti				
L 10mm / Acetal	L 10mm / Titanium				
POI-0021	POI-0024	POI-0031	POI-0033		
Multi-Unit Abutment Analog	Multi-Unit Healing Cap	Multi-Unit Impression-Coping Open Tray	Multi-Unit Impression-Coping Close tray		
L 12mm / Stainless Steel	L 3.5mm / Titanium	L 10mm / Stainless Steel	L 6mm / Stainless Steel		

POI-0046 All our Ball Attachment Angulated (BAA) are supplied with a screw. Tightening the screw to 25Ncm.

POI-0025 All Multi-Unit Accessories are supplied with screw. Tightening the screw to 30Ncm.

# BALL ATTACHMENTS ABUTMENTS



POI-0001	POI-0002	POI-0003	POI-0004	POI-0005		
Ball Attachment H 1mm	Ball Attachment H 2mm	Ball Attachment H 3mm	Ball Attachment H 4mm	Ball Attachment H 5mm		
Titanium	Titanium	Titanium	Titanium	Titanium		
POI-0042	POI-0043	POI-0051	POI-0044	POI-0045	POI-0053	TKI-0033
Ball Attachment 17° / H 2mm	Ball Attachment 17° / H 3mm	Ball Attachment 17° / H 4mm	Ball Attachment 30° / H 2mm	Ball Attachment 30° / H 3mm	Ball Attachment 30° / H 4mm	Key for M.U&balls
Titanium	Titanium	Titanium	Titanium	Titanium	Titanium	Stainless Steel

PATENT PENDING

# BALL ATTACHMENTS ACCESSORIES

POI-0006	POI-0009	POI-0012	POI-0013	POI-0014
Ball Attachment Analog	Metal Housing for Plastic Ball Cap	Soft Plastic Cap Ball	Standard Plastic Cap Ball	Plastic Disc for Ball Attachment
L 12mm / Stainless Steel				


POI-0046 All our Ball Attachment Angulated (BAA) are supplied with a screw. Tightening the screw to 25Ncm.

# EQUATOR



- POI-0056 Equator Attachment H1mm
- POI-0057 Equator Attachment H2mm
- POI-0058 Equator Attachment H3mm
- POI-0059 Equator Attachment H4mm
- POI-0060 Equator Attachment H5mm

 instruments



## TOOLS



### DRILLS

TAG drills are compatible with all TAG implants. The drills are external irrigation, supplied as either standard or short, and drill diameter can be accurately identified through the simple color coding system.



### SURGICAL TOOLS

We produce top quality and attractively designed tools. Kindly note that our tools are not sterilized prior to shipping. The tools must be cleaned, disinfected and sterilized before use. After sterilization, check the packaging of the sterilized instruments for damage. Check the sterilization indicators.

Store tools in a dry, dust-free, vapor-free area at a moderate temperature.

### SURGICAL KITS

TAG offers various types of surgical kits designed to provide complete solutions for implantation procedures.

Kits contain a set of tools and drills that allow implant drilling up to 6.0mm in diameter.





# PREMIUM SURGICAL KIT

KTI-0001



## DRILLS



## SCREW DRIVERS



## PARALLEL PIN



## COUNTER SINK




## RATCHET & DEPTH PROBE



## IMPLANT DRIVERS



## IMPLANT INDICATOR





# COMPACT **axis** SURGICAL KIT









KTI-0006






# COMPACT **massif** SURGICAL KIT









KTI-0002






							
Marking Drill Ø 1.90mm	Pilot Drill External Irrigation Ø 2mm	Pilot Drill External Irrigation Ø 2.5mm	Drill External Irrigation Ø 2.80mm	Conical Drill for implant Ø 3.75mm	Conical Drill for implant Ø 4.2mm	Conical Drill for implant Ø 5mm	Conical Drill for implant Ø 6mm

				
Countersink Ø 3.75mm-Ø 4.2mm	Countersink Ø 5mm-Ø 6mm	Parallel pin L 11mm / S.S	Implant indicator L 14mm / S.S	Ratchet 84mm / S.S

			
Key 1.28mm Long 23mm / S.S	Key 2.42mm Long 23mm / S.S	Motor key 2.42mm Long 28mm / S.S	2.42 Retriever 15mm / S.S, Titanium

							
Marking Drill Ø 1.90mm	Pilot Drill External Irrigation Ø 2mm	Drill External Irrigation Ø 2.80mm	Drill External Irrigation Ø 3.20mm	Pilot Drill External Irrigation Ø 3.80mm	Drill External Irrigation Ø 4.50mm	Drill External Irrigation Ø 5.00mm	Drill External Irrigation Ø 5.50mm

				
Countersink Ø 3.75mm-Ø 4.2mm	Countersink Ø 5mm-Ø 6mm	Parallel Pin Ø 2mm, L 11mm Ø 3.20mm, L 11mm / S.S	Implant indicator L 14mm / S.S	Ratchet 84mm / S.S

			
Key 1.28mm Long 23mm / S.S	Key 2.42mm Long 23mm / S.S	Motor key 2.42mm Long 28mm / S.S	2.42 Retriever 15mm / S.S, Titanium

# PIXEL *crestone* SURGICAL KIT

KTI-0007









-   
 TD1-0002  
 Marking Drill  
 Ø 1.90mm
-   
 TD1-0003  
 Pilot Drill External Irrigation  
 Ø 2mm
-   
 TD1-0004  
 Pilot Drill External Irrigation  
 Ø 2.5mm
-   
 TD1-0005  
 Drill External Irrigation  
 Ø 2.80mm
-   
 TD1-0006  
 Drill External Irrigation  
 Ø 3.20mm






-   
 TK1-0011  
 Square Driver One Piece  
 28mm / S.S
-   
 TK1-0015  
 Ratchet  
 85mm / S.S
-   
 TK1-0012  
 Parallel Pin  
 L 11mm / S.S

# PROSTHETIC KIT

KTI-0008




-   
 TK1-0001  
 Key 1.28mm Short  
 15mm / S. S
-   
 TK1-0002  
 Key 1.28mm Long  
 23mm / S. S
-   
 TK1-0003  
 Key for M.U & Balls  
 S. S
-   
 TK1-0007  
 Motor key 1.28 Short  
 22mm / S. S
-   
 TK1-0008  
 Motor key 1.28 Long  
 28mm / S. S
-   
 TK1-0013  
 Implant indicator  
 L 14mm / S.S

-   
 TK1-0016  
 Torque Ratchet  
 85mm, 15-45Ncm / S.S
-   
 TK1-0006  
 Screw Driver 1.28mm Long  
 25mm / Stainless Steel
-   
 TK1-0030  
 Mock Straight Abutment
-   
 TK1-0031  
 Mock Abutment 15°
-   
 TK1-0032  
 Mock Abutment 25°

## UNIVERSAL HOLDING KIT

KT1-0004




<b>TK1-0034</b> Handle Set 100mm / S. S	<b>TK1-0027</b> Handle Key 1.28mm 43mm / S. S	<b>TK1-0024</b> Extension for 1/4" Holder 50mm / S. S	<b>TK1-0028</b> Handle Key 2.42mm 43mm / S. S
---	---	---	---



## TECHNICIAN HOLDING KIT

KT1-0005



<b>TK1-0034</b> Handle Set 100mm / S. S	<b>TK1-0023</b> 1.28 Kit Screwdriver 140mm / S. S	<b>TK1-0025</b> 1/4" Abutment Fixture Long 15mm / S. S	<b>TK1-0026</b> 1/4" Abutment Fixture Short 13mm / S. S
---	---	--	---




## SCREW REMOVAL KIT

KT1-0009




<b>TK1-0036</b> Screw Removal Burr Guide	<b>TK1-0037</b> 1-72 UNF Tap	<b>TK1-0038</b> Screw Removal Burr
---	---------------------------------	---------------------------------------



\*Drills stoppers available only for standard length  
All drills are External irrigation & made of stainless steel

Diameter	Drill Type	Part Number	Length
2	Short Pilot Drill Ø 2mm	TD1-0011	L 6mm
		TD2-0034	L 6mm
		TD2-0050	L 10mm
	Pilot Drill Ø 2mm	TD1-0003	L 8mm
		TD2-0042	L 8mm
		TD2-0058	L 11.5mm
	Drill Stopper	TD2-0050	L 10mm
		TD2-0058	L 11.5mm
		TD2-0074	L 16mm
2.5	Short Drill Ø2.50mm	TD1-0012	L 6mm
		TD2-0035	L 6mm
		TD2-0051	L 10mm
	Drill Ø2.50mm	TD1-0004	L 8mm
		TD2-0043	L 8mm
		TD2-0059	L 11.5mm
	Drill Stopper	TD2-0051	L 10mm
		TD2-0059	L 11.5mm
		TD2-0075	L 16mm
2.8	Short Drill Ø 2.80mm	TD1-0013	L 6mm
		TD2-0036	L 6mm
		TD2-0052	L 10mm
	Drill Ø 2.80mm	TD1-0005	L 8mm
		TD2-0044	L 8mm
		TD2-0060	L 11.5mm
	Drill Stopper	TD2-0052	L 10mm
		TD2-0060	L 11.5mm
		TD2-0076	L 16mm
3.2	Short Drill Ø 3.20mm	TD1-0014	L 6mm
		TD2-0037	L 6mm
		TD2-0053	L 10mm
	Drill Ø 3.20mm	TD1-0006	L 8mm
		TD2-0045	L 8mm
		TD2-0069	L 13mm
	Drill Stopper	TD2-0053	L 10mm
		TD2-0061	L 11.5mm
		TD2-0077	L 16mm



\*Drills stoppers available only for standard length  
All drills are External irrigation & made of stainless steel

3.8



Drill Length*	TD2-0038	TD2-0046	TD2-0054	TD2-0062	TD2-0070	TD2-0078
	L 6mm	L 8mm	L 10mm	L 11.5mm	L 13mm	L 16mm

4.5



Drill Length*	TD2-0039	TD2-0047	TD2-0055	TD2-0063	TD2-0071	TD2-0079
	L 6mm	L 8mm	L 10mm	L 11.5mm	L 13mm	L 16mm

5



Drill Length*	TD2-0040	TD2-0048	TD2-0056	TD2-0064	TD2-0072	TD2-0080
	L 6mm	L 8mm	L 10mm	L 11.5mm	L 13mm	L 16mm


5.5



Drill Length*	TD2-0041	TD2-0049	TD2-0057	TD2-0065	TD2-0073	TD2-0081
	L 6mm	L 8mm	L 10mm	L 11.5mm	L 13mm	L 16mm




MARKING DRILLS



TD1-0001 Spade Marking Drill Ø 1.50mm  
TD1-0002 Marking Drill Ø 1.90mm

COUNTERSINK DRILLS




TD1-0019 Countersink Ø 3.75 - Ø 4.20mm  
TD1-0020 Countersink Ø 5-Ø 6mm

CONICAL DRILLS




TD1-0021 Conical Drill for implant Ø3.75mm  
TD1-0022 Conical Drill for implant Ø4.20mm  
TD1-0023 Conical Drill for implant Ø5mm  
TD1-0024 Conical Drill for implant Ø6mm



## IMPLANT DRIVERS



## SCREW DRIVERS



# INSTRUMENTS



TK1-0030

Drill Extender  
30mm / Stainless Steel



TK1-0012

Parallel Pin  
Ø 2mm, L 11mm  
Ø 3.20mm, L 11mm/ S.S



TK1-0013

Implant Indicator  
14mm / Stainless Steel



TK1-0015

Ratchet  
85mm / Stainless steel



TK1-0016

Torque Ratchet  
85mm, 15-45Ncm / S.S



TK1-0014

Depth Probe  
Stainless steel



TK1-0020

Abutment Extractor Key  
22mm / Stainless steel



TK1-0021

Open/Close Tray Screw Key  
Stainless steel

# TREPHINE BURRS



TD1-0031

Trepine Burrs  
Ø 3x4mm



TD1-0032

Trepine Burrs  
Ø 4x5mm



TD1-0033

Trepine Burrs  
Ø 5x6mm

TOOLS  
PACKAGE



SUPERSTRUCTURES  
PACKAGE







SyBone®TCP is a porous synthetic ceramic that contains 99% tricalcium phosphate and is designed for the filling of bone voids or defects.

Safety: 100% Synthetic.  
No human or animal origin.

Osteoconductive: Open porosity and highly interconnected with a high mechanical resistance.



0.1mm-0.5mm  
1gr.

ABM-0002  
Granules [0.1 – 0.5mm 1gr]  
indicated for small  
periodontal defects


0.5mm-1mm  
1gr.

ABM-0003  
Granules [0.5 – 1mm 1gr]  
indicated for the filling of cysts  
and the filling of alveolar cavity

1mm-2mm  
1gr.

ABM-0004  
Granules [1 – 2mm 1gr]  
indicated for sinus floor elevation and  
the filling of large cysts





 DIGI TAG



# INNOVATIVE DENTAL SOLUTIONS



## PlanTAG

3D computerized imaging to visualize the anatomy of the patient's jaw. The doctor can plan the exact placement of the virtual implants. This technology provides the doctor with important information about bone quality and quantity.



## GuideTAG

Radiographic template will provide an incredibly high level of accuracy when placing the implants allowing for optimum position, exact depth and angle of the implant.



## ProsTAG

Top down design simplifying the prosthetic stages for immediate loading of customized restorations thanks to prosthetically - driven implant planning and Cad Cam technology.

**tag** dental



**T.A.G Medical Products Corporation Ltd**

Kibbutz Gaaton 2513000 Israel Tel: +972 (0)4 9858400

[www.tag-med.com](http://www.tag-med.com)

---

T.A.G is ISO 13485 and CAN/CSA ISO 13485 certified by LNE  
T.A.G Products are CE marked in accordance with Directive 93/42/EEC

DA-CTLGE Rev.6