# **INSTRUCTIONS FOR USE - DENTAL IMPLANTS**

EN English

FR French

ES Spanish

DE German

RU Russian

UK Ukrainian

IT Italian

HE Hebrew

AR Arabic

It is important to read the Instructions for Use  
in its entirety prior to using the Device

Manufacturer:

**T.A.G. Medical Products Corporation LTD**  
Kibbutz Gaaton 2513000 Israel  
Tel: +972-4-9858400  
Fax: +972-4-9858404

Authorized Representative:

**MEDNET GmbH**

Borkstrasse 10,  
48163 Muenster  
Germany  
Fax: +49 251 3226622

# INSTRUCTIONS FOR USE - DENTAL IMPLANTS

EN

Manufacturer

T.A.G. Medical Products Corporation LTD  
Kibbutz Gaaton 2513000 Israel  
Tel: +972-4-9858400, Fax: +972-4-9858404

EC

Authorized Representative:

MEDNET GmbH

Borkstrasse 10, 48163 Muenster  
Germany  
Fax: +49-251-32266-22

## Indications

The T.A.G. Dental Implant Systems are intended to replace single or multiple teeth in the fully or partially edentulous mandibular or maxillary alveolar process.

The implants are appropriate for immediate loading when good primary stability is achieved and with appropriate occlusal loading.

## Contraindications

TAG Dental implants should not be used in patients who have contraindicating systemic or uncontrolled local diseases such as:

Cardiovascular disorders associated with high endocarditic risk, Coronary insufficiency, Cancers and radiation of the facial region in the past five years, weakened immune system, Diabetes, Bone metabolism disorders, Hemophilia, Poor general state of health, inadequate wound healing capacity, maxillary and mandibular growth not completed, unfavorable anatomic bone conditions, Periodontal disease, Poor oral hygiene, bruxism, Inadequate bone or blood supply, Unrealistic patient expectations or poor patient motivation, Psychological disorders, Smoking, Drug abuse or alcoholism, Steroid use, Allergy to Titanium, Aluminum or Vanadium.

## Precautions

Carefully read these instructions. The sale and use of this device is restricted to, or by the order of, licensed dentist.

## Warnings

- Implants should be placed and restored only by practitioners who are licensed and trained to perform these procedures.
- Adequate preoperative studies should be performed to examine the anatomic structures and to assess the biomechanical, functional, and esthetic requirements of each case.
- Radiographs or other diagnostic reviews should be performed to determine position and topography of the maxillary sinus, nasal cavities, inferior alveolar nerve, mental foramen, natural tooth positions and other anatomical features that may affect implant placement or prognosis.
- Consultation between the surgeon, restorative dentist, and dental laboratory is essential for success.
- Risks of implant placement and restoration include, but are not limited to: infection, implant failure, loss of bone and soft tissue, unfavorable aesthetic result, anesthesia, dysesthesia and paresthesia in the oral and facial areas, sinus infection, dislodgement of implants and instruments in the surrounding structures, damage to adjacent teeth, non-restorable implants, fracture of implants or restorative components, and loosening of implants or restorative components.
- Each implant system has unique measuring characteristics to allow full seating of the implant to the desired depth. In some instances, drill length reference lines measure longer than the stated length of the implant.

It is recommended that the implant surgeon be thoroughly familiar with the specific measurement system being utilized and provide a suitable safety margin adjacent to any teeth and vital structures. Failure to recognize the difference between the actual length of the drill and radiographic measurements can result in permanent injury to the nerves or other vital structures by drilling beyond the depth intended, potentially resulting in permanent numbness to the lower lip and chin or other injuries.

Each implant system has specific design characteristics for mating implants, abutments, prosthetic components, and instrumentation. Combining instruments and components that are not configured or dimensioned for correct mating can lead to mechanical failure of components, damage to tissue, or unsatisfactory aesthetic results.

- The Crestone and 3.3 mm implants are indicated for placement in the incisor region, the Crestone implants are not intended to provide an angle correction.
- One-hundred percent success cannot be guaranteed. Lack of adequate quantity and/or quality of remaining bone, infection, inadequate surgical technique, poor patient oral hygiene, and generalized disease are some potential causes for failure of osseointegration. Both, immediately after surgery or after osseointegration is initially achieved. Pre-operative hard tissue or soft tissue deficits may yield a compromised aesthetic result or unfavorable implant angulation.

With respect to children, routine treatment is not recommended until completion of alveolar growth has been verified.

The TAG Dental implants and accessories have not been evaluated for safety and compatibility in the MR environment. The TAG Dental implants and accessories have not been tested for heating or migration in the MR environment.

All efforts must be made to minimize damage to the host tissue. In particular, special attention must be paid to thermal and surgical trauma and to the elimination of contaminants and sources of infection.

The surgical procedure requires a high degree of precision and care. Any divergence from the principle of least possible trauma at implant installation increases the risk of failure to establish osseointegration.

Please refer to our web site for the specific sequence of drills for each implant type and size. Special care should be taken to avoid over or under preparation of the osteotomy. Implants should be inserted in such a way that they are stable and lack any mobility.

All instruments used in surgery must be maintained in good condition and care must be taken that the instruments do not damage the implants or other components.

Precautions must be taken to avoid the swallowing or aspiration of components used in implant dentistry.

After the implant installation, the surgeon's evaluation of bone quality and initial stability will determine when implants may be loaded.

- An appropriate follow-up protocol should be followed.
- Implants are supplied sterile and are for single use only - DO NOT RESTERILIZE.
- Reprocessing the Implant may result in biocontamination, degraded performance or loss of function.
- Do not use after the expiration date.
- Do not use an implant if it has been opened or the packaging has been damaged. If tampering or damage exists contact your T.A.G. Representative.
- The only allowed modification to the abutments is to shorten the abutment approximately 2mm below to the height of adjacent teeth. Do not shorten the abutment to length less than 4 mm due to the concern of prosthetic retention.
- Plastic cylinder with Ti base are intended to be casted straight only

## Abutment installation procedure

- Remove the healing abutments or cover screws with a hex screwdriver. Now implants are visible and the relevant abutment can be installed.
- Verify that the selected abutment is the one required for installation.
- Insert the abutment into the implant, securing it with a light finger pressure, and fix it with a screwdriver and hex wrench and ratchet torque at 30-35Nm.

## Multi-Unit abutments

The Multi-Unit abutments are intended to provide edentulous patients with an efficient and effective restoration using only four or six implants to support full-arch prosthesis.

Follow casting and burnout instructions as recommended for abutments.  
Multi-Unit plastic cylinder is intended to be casted straight only and connected to the multi-unit abutment.

## Recommended casting procedure

- Waxing Procedure**  
Place the plastic sleeve on top of the analog in the model and secure it with the abutment screw. If necessary shorten the plastic so there is approximately 2mm of clearance with the opposing dentition and approximately 2mm below to the height of adjacent teeth. The minimum abutment length must be 4 mm. Once the correct occlusal height is achieved, proceed with conventional crown and bridge waxing procedures.
- Burnout & casting**  
Follow casting and burnout instructions as recommended by the dental gold alloy manufacturer. The framework should seat passively and be carefully verified. It is recommended that the final restoration be placed using a 30-35Nm torque.

## Follow Up Care

- Patients should be instructed in appropriate oral hygiene and care of the implants and restorations.
- Periodic follow up appointments should be made to confirm and maintain adequate function of the implants and the health of the surrounding tissues.

## CLEANING AND STERILIZATION INSTRUCTIONS

These cleaning and sterilization instructions have been validated for preparing reusable instruments for further use. It is the responsibility of the end user to ensure that the cleaning and sterilization is actually performed using appropriate equipment, materials, and personnel to achieve the desired result. This normally requires validation and routine monitoring of the process. Any deviation from these instructions should be evaluated for effectiveness and potential adverse consequences.

## PREPARATION

- Clean instruments as soon as possible after use. If cleaning must be delayed, wet the instruments in a compatible liquid solution to prevent drying and encrustation of surgical soil.
- Avoid prolonged exposure to saline to minimize the chance of corrosion.
- Remove excessive soil with a disposable wipe.

**Note:** Prior to sterilization, instruments should be cleaned by either automated or manual means as described below.

## MANUAL CLEANING

### Cleaning Instructions

- Prepare an enzymatic cleaning solution in accordance to the manufacturer's instructions.
- Soak soiled instruments for a minimum of 10 minutes. Do not exceed 30°C.
- When cleaning, fully immerse the instrument in the cleaning solution to avoid aerosol generation. Use a soft bristle brush to remove all traces of blood and debris; pay close attention to any hard to reach areas, textured surfaces, or crevices. When cleaning Instruments with cannulations or lumens (i.e. tubes) or holes, use a tight-fitting, soft, non-metallic cleaning brush or pipe cleaner to scrub the cannula, lumen, or hole. Push in and out, using a twisting motion to remove debris.
- Use a syringe filled with enzymatic cleaning solution to flush hard to reach internal areas. When cleaning articulating instruments, (those with moveable parts), brush with a soft non-metallic bristle brush to remove all traces of blood and debris. Pay close attention to threads, crevices, seams, and any hard to reach areas. Actuate any moveable mechanisms, such as hinged joints, box locks, or spring-loaded features, to free trapped blood and debris. If the components of the instrument can be retracted, retract or open the part while cleaning the area.

For instruments with flexible shafts, flex the instrument under the cleaning solution.

### Rinsing Instructions

Rinse the instrument thoroughly with filtered water at 38–49°C, agitate for a minimum of 1 minute; repeat rinse at least two additional times. When rinsing, pay particular attention to flush the cannulations, lumens, or holes with warm tap water. Also pay particular attention to internal areas and moveable parts. Actuate moveable parts while rinsing. If the components of the instrument can be retracted, retract or open the part while rinsing the area.

For instruments with flexible shafts, flex the instrument under the rinse solution.

### Drying Instructions

Dry the instrument immediately after final rinse. Use filtered compressed air to dry internal areas.

## After drying proceed to – INSPECTION AFTER CLEANING INSTRUCTIONS.

## AUTOMATED CLEANING

### General Provisions

- When cleaning, fully immerse the instrument in the cleaning solution to avoid aerosol generation
- When cleaning Instruments with cannulations or lumens (i.e. tubes) or holes, use a tight-fitting, soft, non-metallic cleaning brush or pipe cleaner to scrub the cannula, lumen, or hole. Push

in and out, using a twisting motion to remove debris. Use a syringe filled with enzymatic cleaning solution to flush hard to reach internal areas

- When cleaning articulating instruments, (those with moveable parts), brush with a soft non-metallic bristle brush to remove all traces of blood and debris. Pay close attention to threads, crevices, seams, and any hard to reach areas
- Actuate any moveable mechanisms, such as hinged joints, box locks, or spring-loaded features, to free trapped blood and debris
- If the components of the instrument can be retracted, retract or open the part while cleaning the area. For instruments with flexible shafts, flex the instrument under the cleaning solution
- Process**
- Prepare an enzymatic detergent (ENZOL®) bath – 59.41mL ENZOL®+ 3785mL water (2oz/gallon)
- Make sure that the temperature of the prepared solution is at least 23.5°C
- By using 10mL syringe with a needle, flush each lumen present on the apparatus – for 10 times at least (till no visible debris is present)
- Scrub the cervices, slots and the apparatus with a soft bristled brush until visible debris is fully removed. Pay close attention to any hard to reach areas, textured surfaces, or crevices
- Place the apparatuses in the washer in a manner to facilitate draining
- Make sure that the detergents, which are loaded in washer disinfectant detergent tanks, are suitable for use in automated cleaning cycle (Renuzyme Plus® and Tech Wash III®)
- Use/set the following parameters:

Phase	Time (MM:SS)	Temperature (°C)	Detergent Type
Soaking Stage	00:35	Cold Water	---
Hot Water Wash	02:00	60±5°C	Renuzyme Plus®
Cold Water Wash	00:45	27±5°C	Tech Wash III®
Rinse (cycle 1)	02:00	Hot Water	---
Final Rinse	01:00	82.2±5°C (DI Water)	---
Drying	06:00	82.2±5°C	---

## After the process accomplishment proceed to – INSPECTION AFTER CLEANING INSTRUCTIONS.

## INSPECTION AFTER CLEANING

- Inspect all instruments before sterilization or storage to ensure the complete removal of soil from surfaces, tubes and holes, moveable parts.
- If areas are difficult to inspect visually, check for blood by immersing or flushing the instrument in a 3% hydrogen peroxide solution. If bubbling is observed, blood is present. Rinse instruments thoroughly after using hydrogen peroxide solution.
- If soil is still present, re-clean the instrument.

## MAINTENANCE

Between uses, lubricate moving parts with a water-soluble lubricant in accordance with the manufacturer's instructions.

## INSPECTION AND FUNCTIONAL TESTING

- Visually inspect the instrument and check for damage and wear.
- Cutting edges should be free of nicks and have a continuous edge.
- Jaws and teeth should align properly.
- Moveable parts should have smooth movement without excessive play.
- Locking mechanisms should fasten securely and close easily.
- Long, thin instruments should be free of bending and distortion.

## PACKAGING

- For instruments that are provided in sets, instrument trays may be used for containment during sterilization.
- Package trays/instruments with a barrier wrap material in accordance with local procedures, using standardized wrapping techniques such as those described in ANSI/AAMI ST46.

## STERILIZATION

- The TAG Dental implants are supplied sterile for single patient use.
- Abutment and tools are supplied non-sterile and must be sterilized before use.
- Use a validated, properly maintained and calibrated steam sterilizer.
- Effective steam sterilization can be achieved using the following cycles:

## Cycle Type Temperature Minimum Exposure Time

Pre-vacuum 132–134°C for 4 minutes; Drying time: 15min

Pre-vacuum 134–137°C for 3 minutes; Drying time: 15min

## STORAGE

Store sterile packaged instruments in a manner that provides protection from; dust, moisture, insects, vermin, and extremes of temperature and humidity. For further information related to the use of this instrument, please contact your Representative.

## Materials

TAG implants and abutments are made out of medical grade titanium alloy.

## Disposal

The instruments are to be disposed of as per legislation of the relevant to user country.

## Labeling

Symbols on the product package and labels should be interpreted as follows:

Caution / Attention: See Instructions for Use	Do Not Reuse	Expired by	Manufacturer	Lot Number	Authorized Representative or Distributor Community	Do Not Use If Package Is Open or Damaged	Caution: Federal law restricts this device to sale by or on the order of a physician	Caution: Non-sterile

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MEDNET GmbH, Borkstrasse 10